# Introduction to Multimodal Neuroimaging Structure I: T1, T2, ultra-high field

Howard Aizenstein Maria Ly

Jan 24, 2020





## <u>Outline</u>

Overview of structural MRI research Howard

Review of Key MRI Pulse Sequences Maria

Machine learning Model of Brain Aging
 Maria

Ultra High Field (7T) MR Imaging
 Howard

## MRI Key Concepts





### **Human Subject Issues**

- No radiation exposure (unlike xray, PET, CT),
- thus multiple scans possible


### Risks:


- Loud
- Tight space: concern about claustrophobia
   & participant size
- Metal in/on body can heat, move, or cause artifact

### MR Imaging Key Concepts (Cont)


- Detect signal from <sup>1</sup>H of H<sub>2</sub>O
- In conventional imaging, signal intensity (greyscale) is related to MR relaxation properties of <sup>1</sup>H<sub>2</sub>O such as T1 or T2 relaxation times.
- Water in grey matter has different T1 and T2 relaxation times than water in white matter or CSF. (These are standard structural sequences).
- Some standard sequences (T1, T2, SWI) are discussed today.
- Others include, BOLD fMRI, ASL, DTI, ...







Ultra-High Field (7T) Susceptibility Weighted Imaging






### **MRI Modalities:**

- T1-Weighted ('Anatomic')
  - Morphometry: regional volumetrics, cortical thickness, shape
- T2-Weighted (with or without FLAIR)
  - WMH, lacunes
- SWI (Susceptibility weighted Imaging)
  - Microbleeds, perivascular spaces, vessels, anatomic landmarks

## T1-Weighted Anatomic Scans

- Weighted towards imaging differences in T1 (longitudinal) relaxation properties of tissue.
- Short TR, short TE,
- Accentuates gray/white contrast
- Often Volumetric High resolution, scan such as SPGR, MPRAGE is used



1.5 Tesla 3 Tesla 7 Tesla

All rely to some extent on basic image processing functions of


**Registration** (pair-wise): align one image to another


Segmentation: label voxels which meet a criteria

- Atlas-Based Segmentation
- Voxel-Based Morphometry (VBM)
- Tensor-Based Morphometry (TBM)
- Cortical Thickness
- Shape Analyses

- Atlas-Based Segmentation
  - Registering the atlas brain into the individual brain space.
  - Specific regions in the individual brains are then segmented using the computed spatial transformation and the MNI anatomical template.

### **Cross Registration**





Atlas-Based Segmentation:

template

- Voxel-Based Morphometry (VBM)
- Tensor-Based Morphometry (TBM)
- Cortical Thickness
- Shape Analyses

All rely to some extent on basic image processing functions of **Registration** (pair-wise): align one image to another **Segmentation**: label voxels which meet a criteria

Atlas-Based Segmentation:

template

- Voxel-Based Morphometry (VBM):
- Tensor-Based Morphometry (TBM):
- Cortical Thickness
- Shape Analyses

All rely to some extent on basic image processing functions of

Registration (pair-wise): align one image to another

**Segmentation**: label voxels which meet a criteria




Atlas-Based Segmentation:



- Voxel-Based Morphometry (VBM):
- Tensor-Based Morphometry (TBM):

temmate




Atlas-Based Segmentation:

template individual

Voxel-Based Morphometry (VBM):

• Tensor-Based Morphometry (TBM):

template



Atlas-Based Segmentation:

template

- Voxel-Based Morphometry (VBM):
- Tensor-Based Morphometry (TBM):
- Cortical Thickness
- Shape Analyses



• Atlas-Based Segmentation:

template




- Voxel-Based Morphometry (VBM):
- Tensor-Based Morphometry (TBM):
- Cortical Thickness

individual



• Shape Analyses




Atlas-Based Segmentation:

template

individual

temMate


- Voxel-Based Morphometry (VBM):
- Tensor-Based Morphometry (TBM):
- Cortical Thickness
- Shape Analyses




Tamburo, Aizenstein, et al., 2007, ISBI

## T2-Weighted (with or without FLAIR)


- Weighted towards imaging differences in T2 (transverse) relaxation properties of tissue.
- Long TR, short TE
- Often used with FLAIR to accentuate WMH/CSF contrast



T1-Weighted




T2-Weighted



T2-Weighted FLAIR

#### White Matter Hyperintensities (WMH)

- Areas of hyperintensities found on MRI images, best seen in T2-weighted FLAIR images.
- Located around ventricles (periventricular) and in deep WM.
- Caused by white matter alterations such as demyelination, gliosis, or axonal degeneration
- Most studies point to a vascular (ischemic) origin for WM changes



#### Clinical Significance

- Some studies have shown an association between WMH and disorders of aging including
  - Cognitive deficits
  - Depression

## Vascular Depression Hypothesis of LLD

Cerebral small vessel disease contributes to the occurrence and perpetuation of depression in some older adults. (Alexopoulos, 1989)

#### **Evidence:**

- WMHs (subcortical) associated with LLD (Hickie et al, 1997; Krishnan, 1997)
- DTI shows WM changes assoc. with poor response. (Alexopoulos, et al., 2002)
- Decreased resting CBF/metabolism in PFC, caudal ACC, striatum in LLD (Nobler et al., 1999; Meltzer et al. 1998)
- LLD show impairment on executive and speed tasks associated with poor response (Alexopoulos et al., 2000)

### Measurement of WMH

Two approaches used to evaluate WMH on FLAIR images:


### Visual rating systems

- WMH severity is visually graded by trained experts
- · High variability based on rater and scales

### Computer-assisted volumetric evaluation


- More reliable than visual ratings
- Allows for accurate localization and volumetric analysis of WMH

### **Automated WMH Segmentation**



- Fuzzy connected algorithm to segment the WMH
- Left-side is unprocessed FLAIR, right-side is segmented WMH


## Distribution of WMH with Depressive Symptoms



From Health ABC study, Rosano, Aizenstein, Tudorascu

## Strokes (including Lacunes)

Some partially automated methods, generally done manually.



From: Silent Brain Infarcts: A Review of MRI Diagnostic Criteria; Zhu et. al. http://stroke.ahajournals.org/content/42/4/1140.full

## Population Neuroscience

How is MRI research relevant for population studies?

Isn't it only for small sample, mechanistic studies?


- Cost isn't prohibitive.
- Number of population based studies with 100's and 1000's of participants. CHS, Health ABC, WHI, BLSA, ADNI

How can it guide public health?

- Examples: New biomarkers for AD diagnosis, concept of preclinical AD.
- MRI can serve as proxy-outcome for intervention studies.
- And as a neuroendophenotype for behavioral disorders

## Diffusion Weighted Imaging

### Diffusion or Brownian Motion




Hagmann, Radiographics, 2006


- Attenuation of the signal due to diffusion is acquired using pulsed excitation gradient in different directions.
- Used clinically as early indicator of stroke (cell swelling)
- Advanced methods such as DTI, DSI, HDFT used to assess myelin integrity and do fiber tracking.

## Diffusion Imaging

• Raw diffusion images for detecting stroke



## Diffusion Tensor Imaging (An)isotropic diffusion



Mukherjee et.al., AJNR (2008) 29, 632-641

Fractional Anisotropy: a measure of directional diffusivity.

Higher □ greater myelin integrity.




## Highlights of the Translational Studies in the 7 Tesla Bioengineering Research Program

All of the 7 Tesla data were obtained using the Tic Tac Toe RF Coil System (conceived and developed in the RF Research Facility).

Tamer S. Ibrahim (PI)

http://rf-research-facility.engineering.pitt.edu





### MRI Modalities:

- T1-Weighted ('Anatomic')
  - Morphometry: regional volumetrics, cortical thickness, shape
  - Perivascular spaces
- T2-Weighted (with or without FLAIR)
  - WMH, lacunes
  - Anatomic landmarks
- SWI (Susceptibility weighted Imaging)
  - Microbleeds, perivascular spaces, vessels, anatomic landmarks
- TOF, CSF pulsatility

## UHF (7T) MRI

- Increased SNR due to aligning more protons
- Challenge Larmor wavelength through tissue at 7T is 11 cm
- Leads to inhomogeneities.
- RF solution with Tic-Tac Toe design






NIH R01 MH111265 NIMH and R01 EB009848 NIBIB http://rf-research-facility.engineering.pitt.edu

### T1-weighted **MPRAGE**

- Very crisp well-defined gray-white differentiation is achieved at 7T. Able to resolve/visualize small perivascular spaces better at 7T than 3T.
  Radiofrequency field 7T inhomogeneities are mitigated using the Tic Tac Toe RF coil system with extended RF field coverage to temporal lobe & cerebellum


NIH RF1/R01 AG025516 NIA

**Normal Aging Study (3T)** 




**Normal Aging Study (7T)** 

T1 MPRAGE 0.55mm x 0.55mm x 0.55mm TA = ~10min TR = 3.65s; TE = 2.53ms; TI = 2.9s



**T1 MPRAGE** 1mm x 1mm x 1.2mm TA = ~9minTR = 2.3s; TE = 2.98ms; TI = 2.5s


### Enhanced GM contrast at 7T



RF Research Facility

http://rf-research-facility.engineering.pitt.edu

•Well visualized substructure anatomy of the hippocampus at 7T.
•Radiofrequency (RF) field inhomogeneities are mitigated using the 7T Tic Tac Toe RF coil system.
•Extended RF field coverage to temporal lobe & cerebellum even W. RF Power-Intensive T2 scan.
Small Vessel Disease/Late-Life Depression Study (7T) Small Vessel Disease/Late-L Small Vessel Disease/Late-Life Depression Study (3T)



2D TSE 0.375mm x 0.375mm x 1.5mm TA = ~3.5min TR = 10.06s; TE = 61ms

http://rf-research-facility.engineering.pitt.edu

**3D SPACE** 0.8mm x 0.8mm x 0.8mm


TA = ~6min

TR = 3.2s; TE = 563ms

## 7T Allows for Accurate and Highly Detailed Automatic Hippocampus


http://rf-research-facility.engineering.pitt.edu


**Segmentation** (Midlife Neurocognitive Disparities: NIH R01 DK110041 NIDDK)



- Combination of phase and magnitude, enhances 'susceptibility signal'
- Contrast scales very well with increasing field strength
- Identifies microbleeds: prominent in cerebral amyloid angiopathy and increasingly recognized as important marker of cerebrovascular integrity in elderly.
- Images venous side of vascular tree.
- Quantitative Iron.


Contrast scales very well with increasing field strength





7T SWI (collaboration with Drs. Rosano and Boada)

 Identifies microbleeds: prominent in cerebral amyloid angiopathy and increasingly recognized as important marker of cerebrovascular integrity in elderly.



**7T SWI** 


Images venous side of vascular tree (deoxygenated hemoglobin)

RF Research Facility http://rf-research-facility.engineering.pitt.edu

Zoomed previous slide: Much better visualization (at 7T) of medullary veins in the deep periventricular white matter (see many more vessels, also the vessels are better defined).
 White matter lesion with a traversing vein (can not be visualized at 3T). White matter lesion without any convincing traversing vein (can not be visualized at 3T). There is value in discerning these two subtypes of white matter lesion to distinguish traditional small vessel disease from an inflammatory component.

NIH RF1/R01 AG025516 NIA **Normal Aging Study (3T)** 

**Healthy Brain Initiative Study (7T)** 




0.22mm x 0.22mm x 1.5mm TA = ~10min TR = 1.96s; TE = 14ms

2D SWI

0.5mm x 0.5mm x 1.5mm TA = ~9minTR = 28ms; TE = 20ms


3D SWI

# Susceptibility Weighted Imaging (0.2 x 0.2 x 1.5 mm³) shows cortical microvessels at sub 100 microns *in-vivo*



Voxel containing long

medullary arterioles








e or superimposed capillal 4.) Bar = 500 μm.




Fig. 2.—Arteriolar supply to cortex Light photomicrograph of 100-μm celloidin section with alkaline phosphates vascular stain contertsained with cresty violet acetate and light green. Medium-stzed pial artery and vein are encased in common connective tissue sheath in space between two gyri. Cortical arterioles penetrate cerebral surface. It is possible to differentiate gray matter from white matter by density of capillary bed. This is an example of single supply with adjacent arterioles arising from a common source. Cortical arterioles share capillary network. Note capillary-free zone (arrows) surrounding alkaline phosphatase-stained arterioles, which is anatomic evidence that arterioles have a primary nutritive responsibility to their surrounding neural tissues. This was observed [4] prior to general recognition of this fact by physiologists [5]. This periarteriolar capillary-free zone will not be easily distinguished on thicker specimens imaged with microradiography because of superimposed capillaries above and below particular arteriole. (See Fig. 4.) Bar = 500 μm.

### **BrainVein**

- Computer Vision algorithm for Segmentation of medullary small veins from UHF SWI.
- Takes advantage of known radial orientation of vessels



Segmented 3D veins



Stetten et al., CARS 2019


G. Stetten, Minjie Wu.


#### Role of CSF Clearance in late-onset AD

- Amyloid and other neurotoxic metabolites are cleared through CSF
- Decreased clearance is salient model for late-onset AD
- The CSF is formed from interstitial fluid pushed along small channels from pulsations in the vessels
- Related to delta-waves (slow-wave sleep).
- Team of investigators addressing this including bioengineers and sleep researchers.

**BIOE RF Lab:** Tamer Ibrahim. Tales Santini, Tiago Martins

Kristine Wilckens, Minjie Wu, Caterina Rosano, Joe Mettenburg








# 100ms TR: Bold Contrast


**BIOE RF Lab:** Tamer Ibrahim. Tales Santini, Tiago Martins

# And Pulsatility




**BIOE RF Lab:** Tamer Ibrahim. Tales Santini, Tiago Martins


### High resolution MR angiography without contrast



**BIOE RF Lab:** Tamer Ibrahim. Tales Santini, Tiago Martins

RF Research Facility http://rf-research-facility.engineering.pitt.edu





 $(0.375x0.375x0.375 \text{ mm}^3)$  Whole Head (15-cm), Total scan time ~ 12min

### **Future Directions**

A current research challenge involves relating the different image modalities, essential for translational research.

#### Mechanism:

Can UHF imaging identify the cellular or microstuctural features that explain altered function?

#### Biomarkers:

Statistical approaches, including machine learning, combine across different modalities to better identify predictors of treatment response. (personalized medicine)